

DirectVision Equipment Checklist

Standard Equipment for Performing a DirectVision Catheterization

- DirectVision System (including DirectVision Camera Hand-piece)
- DirectVision Visual Guide (sterile and wrapped from SPD)
- DirectVision Catheterization Kit
- DirectVision Catheters (14fr/18fr)
- Sterile NS for Irrigation
- Lidocaine Injection (UroJet) if applicable
- Extra set of Sterile Glove

3rd Party Equipment and Accessories for Advanced Use: CBI or Guidewire

- Cook Medical Foley Side-Arm Adapter (#050202) for Guidewire insertion
- Hydrophilic Guidewire such as a Boston Scientific Zipwire® (Straight #630-205B 15CM by .035in)
- Standard Stopcock Valve or Female Foley Straight Adapter for CBI
- StatLock 3-Way Foley for securing drainage bag (Bard Medical #FOL0105)

R002